









PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE. Always swab self-sealing port with sterile alcohol prior to accessing with a sterile syringe. The NR Protocol requires an additional 20mL syringe that is not included in the kit.

<p>Step 1:</p>  <p>Using the filter needle, draw up 2mL of Citrate Anticoagulant into the 20mL Syringe.</p>	<p>Step 2:</p>  <p>Draw 18mL whole blood from the patient, filling the syringe to 20mL.</p>	<p>Step 3:</p>  <p>Remove the red cap and slowly load the 20mL of anticoagulated whole blood into the Concentrating Device. Recap the device using the clear cap (make sure to remove the cap cover).</p>
<p>Step 4:</p>  <ol style="list-style-type: none"> Place the Concentrating Device and Counterbalance Device* in the centrifuge buckets at opposite ends. Process at: <div style="border: 1px solid black; padding: 5px; margin: 10px 0;"> <p>Executive Series Centrifuge 5 minutes - 4400 RPM</p> <p>Platinum Series Centrifuge AbsolutePRP/BMC</p> </div> <p>*Counterbalance Device needs to be filled with water weighing equal to the Concentrating Device (+/- 3g).</p>	<p>Step 5:</p>  <p>Piston touches the RBC interface (trapping plasma & platelet buffycoat inside)</p> <p>Attach the 20mL syringe* and aspirate the plasma until the piston touches the RBC interface (trapping the plasma & platelet buffycoat inside) then stop aspirating. Remove the 20mL syringe.</p> <p>*Extra 20mL syringe for LR protocol is not included with kit</p>	<p>Step 6:</p>  <p>3-4mL PRP (Plasma & platelet buffycoat)</p> <p>Attach the 12mL syringe and aspirate 3-4mL of PRP (Gently agitate the syringe to re-suspend the platelet buffycoat into the plasma).</p>

If you have any questions, please reach out to our Clinical Education Team:
education@plymouthmedical.com or 888.392.5076 Ext. 3