GS120 PUREPRP TWO NP-PRP IFU PLYMOUTH MEDICAL



PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE. Always swab self-sealing port with sterile alcohol prior to accessing with a sterile syringe.






Using the filtered needle, draw up 6mL of Sodium Citrate Anticoagulant into each 60mL syringes.

Step 2:



Using the butterfly needle draw 54mL whole blood in each syringe from the patient, filling each syringe to 60mL. Collect 120mL of solution.

Step 3:



Load anticoagulated whole blood into each Separator Device.

Step 4:




Make sure each device has equal volume and place at opposite ends in the centrifuge rotor.

> Sapphire Series Centrifuge: Set to PUREPRP 60 SPIN 1

> Platinum Series Centrifuge: Set to PUREPRP SP SPIN 1


> **Executive Series Centrifuge:** Set to 2 minutes and 4400

Step 5: For each Separator Device do the following:



Connect the 60mL syringe to the top port and aspirate the platelet plasma suspension (PPS). Optionally aspirate additional 0.5mL of RBC for optimal platelet recovery.

Step 6:





PLACE CLEAR CAP ON BOTTOM PORT

Transfer all of the platelet plasma suspension (PPS) collected from both Separator Devices into the single Concentrator Device.

> PLACE CLEAR CAP ON BOTTOM PORT PRIOR TO CENTRIFUGATION.

WEB - PLYMOUTHMEDICAL.COM

GS120 PUREPRP TWO NP-PRP IFU PLYMOUTH MEDICAL









IMPORTANT! MAKE SURE CLEAR CAP IS PLACED ON BOTTOM PORT PRIOR TO CENTRIFUGATION!




Counterbalance with equal volume then place directly opposite to each other in the centrifuge rotor buckets.

> Sapphire Series Centrifuge: Set to PUREPRP 60 SPIN 2

> Platinum Series Centrifuge: Set to PUREPRP SP SPIN 2

Executive Series Centrifuge: Set to 4 minutes and 4400


Step 8:



After centrifugation:

Platelet concentrate buffycoat separates out at the bottom of the Concentrator Device.

Step 9:



Attach 60mL syringe to bottom port. Aspirate platelet poor plasma from the Concentrator Device.

Leave 14mL or the desired amount.

Step 10:




Attach the 20mL syringe to bottom port and gently swirl to resuspend the platelet buffycoat into the plasma.

Step 11:



Tilt to immerse the Aspirating Pipe into the PRP.

Step 12:



Extract the PRP into the 20mL syringe.

If you have any questions, please reach out to our Clinical Education Team: education@plymouthmedical.com or 888.392.5076 Ext. 3