



NITRONOX PLUS[®] 0-70%

FDA Cleared Nitrous Oxide and
Oxygen Analgesia System

Help your patients relax

Nitrous oxide is a safe and effective sedative agent that is mixed with oxygen. Nitrous oxide is used during procedures and treatments to help your patients feel more comfortable.

HOW DOES THE NITRONOX PLUS[®] WORK?

The Nitronox Plus[®] 0-70% system delivers an adjustable N₂O/O₂ mixture (0-70%) which is inhaled through a mouthpiece or mask, providing your patients control through self-administration. The patients breathe normally through their mouth, and within a few short minutes they will start to feel the relaxing, anxiolytic effects of the gas mixture. The effects of nitrous oxide wear off within 5 minutes after the last inhalation of the N₂O/O₂ mixture.

APPLICATION RECOMMENDATIONS

The Nitronox Plus[®] 0-70% system is ideal for short, minimally invasive and non-invasive procedures where there is a need for pain management or minimal sedation:

- Bone Marrow Harvesting
- Spine Procedures
- Liposuction
- IV Starts/Blood Draws
- Botulinum Injections
- Childbirth
- Joint Injections
- Incision and Drainage
- Sutures
- Lumbar Punctures
- Cosmetic Procedures
- CT Scans
- Vein Procedures
- Fracture Reductions
- Intraosseous Injections
- Hair Injections & Grafting
- Urogenital Procedures

BENEFITS

Rapid Onset & Short-Acting

- Demand flow option is safe and provides patients with control over pain mitigation
- Reduce recovery time
- Improve throughput

Improved Patient Care

- Relaxes and provide comfort for your patients
- Improves patient satisfaction scores
- Get patients back to normal activity faster

Cost-Effective

- Low consumable cost
- Reduce staffing for patient monitoring
- Complete procedures more efficiently





COMPETITIVE ADVANTAGES

1. The Nitronox Plus® 0-70% is made in the USA and built to last.
2. The Nitronox Plus® 0-70% has never had a recall; 510K FDA cleared K202480.
3. The Nitronox Plus® 0-70% is the only available N₂O & O₂ demand flow system with adjustable N₂O.
4. The Nitronox Plus® 0-70% device can be set to provide 100% oxygen to patients to flush their system of any residual nitrous oxide at the end of the procedure to help patients recover faster.
5. The Nitronox Plus® 0-70% offers additional safety features compared to similar systems on the market (Emergency Air Intake, Key Lock, Function Test Button).

TECHNICAL + SAFETY INFORMATION

- **Oxygen Piloted Regulator (Failsafe)** – prevents flow of N₂O without O₂
- **Dual Seal Diaphragm Block** – reduces risk for leak or failure
- **Patient Self-Administered** – provides control over administration as the patient handheld supply valve is only activated by an inspiratory effort (inhalation); patient can stop inhaling N₂O /O₂ mixture at any point
- **Visual Pressure Gauges** – confirm gas supply with three visual pressure gauges - N₂O, O₂ + gas mixture
- **Audible Alarm** – alerts staff when either O₂ or N₂O gas supply pressure drops to approx. 35 psi (2.41 bar)
- **Nitrous Oxide Regulator Check Valve (for 2-cylinder cart)** – prevents flow back out into the room from the hose when removing N₂O regulator
- **Oxygen Enrichment Feature** – automatically increases O₂ concentration by an additional 15% in case of shallow breathing
- **Emergency Air Intake** – allows patient to pull in room air through breathing circuit should gas supply run out
- **100% Oxygen Capability** – flushes their system of any residual nitrous oxide at the end of the procedure
- **Internal, Accessible Demand Valve** – offers security and quick access for in-field servicing
- **Indexed Gas Line Connections and Cylinder Regulators** – prevents cross connection
- **Key Lock System** – renders device unusable and helps to deactivate device for compliance and safety
- **Function Test Button** – provides easy access for periodic safety and function tests
- **Optional Built-In Patented Scavenging Interface** – ensures the safe, effective removal of exhaled nitrous oxide, if needed
- **Quantifiable Exposure** – ensures certainty over N₂O exposure with our available monitoring badges

Regulatory: K202480
Class II breathing gas mixer

TEL – 888.392.5076
FAX – 888.392.1455
WEB – PLYMOUTHMEDICAL.COM

PLYMOUTH MEDICAL
135 PLYMOUTH ST, SUITE 410
BROOKLYN, 11201 NY