

Nitronox Plus® 0-70% Key Features



PLYMOUTH
MEDICAL

Emergency air intake valve – allows patient to pull in room air through breathing circuit should gas supply run out

3 visual gas pressure gauges – real time confirmation of gas supply and gas mixture pressure

Built-in oxygen failsafe – prevents the flow of nitrous oxide without oxygen



Mixed gas outlet – breathing circuit connects here

Single adjustable nitrous oxide knob – easily control the percentage of nitrous oxide (adjustable models only)

On/off key lock for security – nitrous oxide will only flow in the on position, oxygen flows in either position

Audible low gas pressure alarm – alert indicates either low oxygen or nitrous oxide gas pressure

Diameter indexed safety system (DISS) connections – prevents crossing of gas lines



Demand valve access panel – offers quick access for in-field servicing

Built in demand valve – designed for both pediatric and adult patients

Function test button – easy access to perform periodic safety and function tests