



GS30 PurePRP II (30mL Double Spin Kit)	PRP Volume	PLT Concentration	PLT Dose*	PLT/mL
	1mL	28x	7,000,000,000	7,000,000,000
	2mL	14x	7,000,000,000	3,500,000,000
	3mL	9.3x	7,000,000,000	2,333,333,333
<b>Suggested Protocol</b>	4mL	7.0x	7,000,000,000	1,750,000,000
	5mL	5.6x	7,000,000,000	1,400,000,000
	6mL	4.7x	7,000,000,000	1,166,666,667
	7mL	4.0x	7,000,000,000	1,000,000,000
	8mL	3.5x	7,000,000,000	875,000,000

GS60 PurePRP II (60mL Double Spin Kit)	PRP Volume	PLT Concentration	PLT Dose*	PLT/mL
	1mL	49x	12,250,000,000	12,250,000,000
	2mL	24.5x	12,250,000,000	6,125,000,000
	3mL	16.3x	12,250,000,000	4,083,333,333
	4mL	12.3x	12,250,000,000	3,062,500,000
	5mL	9.8x	12,250,000,000	2,450,000,000
	6mL	8.1x	12,250,000,000	2,041,666,667
<b>Suggested Protocol</b>	7mL	7.0x	12,250,000,000	1,750,000,000
	8mL	6.1x	12,250,000,000	1,531,250,000
	9mL	5.4x	12,250,000,000	1,361,111,111
	10mL	4.9x	12,250,000,000	1,225,000,000
	11mL	4.5x	12,250,000,000	1,113,636,364
	12mL	4.1x	12,250,000,000	1,020,833,333
	13mL	3.8x	12,250,000,000	942,307,692
	14mL	3.5x	12,250,000,000	875,000,000

GS120 PurePRP II (120mL Double Spin Kit)	PRP Volume	PLT Concentration	PLT Dose*	PLT/mL
	1mL	98x	24,500,000,000	24,500,000,000
	2mL	49x	24,500,000,000	12,250,000,000
	3mL	32.7x	24,500,000,000	8,166,666,667
	4mL	24.5x	24,500,000,000	6,125,000,000
	5mL	19.6x	24,500,000,000	4,900,000,000
	6mL	16.3x	24,500,000,000	4,083,333,333
	7mL	14x	24,500,000,000	3,500,000,000
	8mL	12.3x	24,500,000,000	3,062,500,000
	9mL	10.8x	24,500,000,000	2,722,222,222
	10mL	9.8x	24,500,000,000	2,450,000,000
	11mL	8.9x	24,500,000,000	2,227,272,727
	12mL	8.2x	24,500,000,000	2,041,666,667
	13mL	7.5x	24,500,000,000	1,884,615,385
<b>Suggested Protocol</b>	14mL	7x	24,500,000,000	1,750,000,000
	15mL	6.5x	24,500,000,000	1,633,333,333
	16mL	6.1x	24,500,000,000	1,531,250,000
	17mL	5.8x	24,500,000,000	1,441,176,471
	18mL	5.4x	24,500,000,000	1,361,111,111

\*Assuming patient baseline of 250,000 PLT/ uL