



GS301-EU Device Comparison



PLYMOUTH
MEDICAL

	Existing Device AbsolutePRP™ Gold 20mL	Upgraded Device AbsolutePRP™ Gold 20mL	Comparison Outcome
			
Final Output Volume	7-11mL	7-11mL	Unchanged
Single Spin Processing	Yes	Yes	Unchanged
Processing Settings	Executive Series 4000 RPM / 1 min Platinum Series AbsolutePRP Gold	Executive Series 4000 RPM / 1.5 min Platinum Series PUREPRP SPIN 1	Upgraded
Self-Sealing Aspirating Port	No	Yes	Upgraded
Aspirating Pipe	Yes	No	Upgraded
Closed System Processing	No	Yes	Upgraded
Anticoagulant	Anticoagulant ACD-A	Anticoagulant ACD-A	Unchanged
> 90% Platelet Yield	Yes	Yes	Unchanged
2X to 3X Platelet Concentration	Yes	Yes	Unchanged
Centrifuge Bucket Insert Required	Yes	No	Upgraded

If you have any questions, please reach out to our Clinical Education Team: education@plymouthmedical.com or 888.392.5076 Ext. 3

TEL - 888.392.5076
FAX - 888.392.1455
WEB - PLYMOUTHMEDICAL.COM

PLYMOUTH MEDICAL
135 PLYMOUTH ST, SUITE 410
BROOKLYN, NY 11201

V3-11-22