

# NucleoCounter® - NC-200™ Installation Checklist



PLYMOUTH  
MEDICAL

Please check each box to indicate the item is available and ready for installation. Return this form via e-mail at least 3 business days before the scheduled installation date. If any item is not available, please inform us or your Chemometec field application scientist ASAP.

## The NucleoCounter® instrument includes:

- Laptop stand (fits a 14 inch laptop)
- NucleoCounter® instrument power supply and cable
- Chemometec USB flash drive (the flash drive is under the yellow lid)
- USB cable

## Required items which are not included with instrument purchase:

- Desktop computer or laptop with administrative login information or Windows tablet (If it's a Windows 7 PC, please have latest Windows Updates installed)
- Via1-Cassettes - available through Plymouth Medical [here](#)
- Solution 17 - Blood Lysis Buffer 25 ml - available through Plymouth Medical [here](#)
- Pipette 20-200µL - recommendation [here](#)
- Pipette Tips 20-200µL - recommendation [here](#)
- Eppendorf Microtubes 1.5mL - recommendation [here](#)
- Incubator for Eppendorf Microtubes 1.5mL - recommendation [here](#)

There may be additional items that are required based on the type of quantification you would like to perform with the NucleoCounter®. We will let you know at the time of the instrument purchase if you require additional items.

## Optional items:

- Computer Mouse

If you have any further questions, please don't hesitate to contact us via email [customerservice@plymouthmedical.com](mailto:customerservice@plymouthmedical.com) or phone (888) 392-5076 Ext. 3.